



CARDIOVASCULAR  
RESEARCH  
TECHNOLOGIES

# Zilver PTX Drug-Eluting Stent Mortality Analysis

Gary Ansel, MD, FACC  
OhioHealth Heart and Vascular Physicians  
Columbus, OH

***iIMPACT*** YOUR PRACTICE

## Gary Ansel, MD, FACC

### Disclosures:

Advisory and consultation services: Medtronic, CR Bard, Boston Scientific, Cook Medical, Philips Medical, WL Gore, Veryan Medical, Reflow Medical, Shockwave, Intact Vascular, Abbott Vascular, Alucent, Contego Medical, Cardiovascular Systems. Surmodics, Vascular Dynamics, Vatrix, and VIVA Physicians.


# Recent Correction to 5-year Zilver PTX Publication

- Katsanos K, et al. meta-analysis published December 6, 2018 in JAHA
- Data reviewed and errors identified in 5-year Zilver PTX publication
  - Incorrect patient flow diagram submitted during final publication process
  - Mortality numbers transposed in overall primary randomization comparison
- Corrections submitted to Circulation on December 18, 2018 and published on February 19, 2019

| Risk Ratio (95% CI) for<br>All-cause death at 4 to 5 years |                     |
|--|---------------------|
| Based on original figure                                   | 1.94 (1.28 – 2.96)* |
| Based on corrected figure                                  | 1.66 (1.14 – 2.44)  |


\* Katsanos K, et al. 2018. JAHA

# RCT Patient Flowchart



# Zilver PTX RCT

## 5-year Mortality Analysis



PTA\*  
n = 237  
Died = 24  
KM = 15.3%

Zilver PTX  
n = 242  
Died = 41  
KM = 22.1%


**p=0.04**

**\*40% of PTA group = Zilver PTX**

# Zilver PTX Key Points


- Data available to Katsanos K, et al. did not identify all patients who were treated with a Zilver PTX stent
  - Patient-level data were not used in the analysis
  - 40% of patients in the PTA group were treated with a Zilver PTX stent
- Patient level analysis demonstrates no difference in mortality rate for Zilver PTX compared to PTA/BMS
  - Causes of death for Zilver PTX are similar to PTA/BMS

# RCT Patient Flowchart




The RCT study design allowed optimal PTA patients requiring reintervention within the first year post-procedure to cross over to treatment with the Zilver PTX stent

# RCT Patient Flowchart






# PTA Group Composed of Zilver PTX Patients



**40% of PTA group = Zilver PTX**  
**70% of patients in study = Zilver PTX**

# Zilver PTX RCT

## 5-year Mortality Analysis



PTA / BMS

n = 143

Died = 17

KM = 17.6%

Zilver PTX

n = 242

Died = 41

KM = 22.1%

?


Zilver PTX

n = 94

Died = 7

KM = 9.4%

# Zilver PTX RCT Final 5-year Mortality Analysis



PTA / BMS  
n = 143  
Died = 17  
KM = 17.6%

Zilver PTX  
n = 336  
Died = 48  
KM = 18.7%

**p=0.53**

**No significant difference  
between Zilver PTX  
and PTA / BMS**

# Covariate Analysis – RCT





- Cox proportional hazards model
- Included comorbidities that may be related to mortality as well as other factors of interest
- No significant difference between Zilver PTX and PTA / BMS ( $p=0.51$ )

| Covariate                 | Multivariate p-value |
|---------------------------|----------------------|
| Age                       | <b>0.0002</b>        |
| Congestive heart failure  | 0.09                 |
| Diabetes                  | 0.11                 |
| Lesion length             | 0.12                 |
| Carotid disease           | 0.13                 |
| Claudication/CLI          | 0.14                 |
| Smoking                   | 0.17                 |
| Cardiac arrhythmia        | 0.21                 |
| Hypertension              | 0.46                 |
| Gender                    | 0.50                 |
| PTX vs. PTA/BMS           | 0.51                 |
| Country (US, JP, Germany) | 0.59                 |
| Pulmonary disease         | 0.61                 |
| Hypercholesterolemia      | 0.63                 |
| Previous MI               | 0.99                 |

# Dose Analysis

- Meta-analysis from Katsanos incorrectly identified Zilver PTX as a high dose device
  - Total amount of paclitaxel on a Zilver PTX stent is approximately 10% to 20% of the amount on a DCB
- Zilver PTX has similar total amount of paclitaxel compared to Eluvia with no polymer and a shorter paclitaxel exposure

# Dose Analysis

| Device                      | Paclitaxel Density                              | Total Paclitaxel Load<br>(7 x 80 mm) |   | Paclitaxel Exposure                |
|-----------------------------|---|--------------------------------------|---|------------------------------------|
| Boston Scientific<br>Eluvia | 0.167 $\mu\text{g}/\text{mm}^2$<br>total area   | 0.3 mg                               |  | $\geq 1$ year<br>permanent polymer |
| Cook<br>Zilver PTX          | 3 $\mu\text{g}/\text{mm}^2$<br>abluminal area   | 0.7 mg                               |  | 2 months<br>polymer free           |
| Bard<br>Lutonix DCB         | 2 $\mu\text{g}/\text{mm}^2$<br>abluminal area   | 3.5 mg                               |  | < 2 months                         |
| Medtronic<br>In.Pact DCB    | 3.5 $\mu\text{g}/\text{mm}^2$<br>abluminal area | 6.9 mg                               |  | < 2 months                         |

References: Device SSEDs/IFUs; Müller-Hülsbeck, Expert Opinion on Drug Delivery 2016, Dake, et al. JVIR 2011; Gongora, et al. JACC Cardio Interv, 2015; <http://www.bostonscientific.com/en-US/products/stents--vascular/eluvia-drug-eluting-stent-system/sustained-drug-release.html> (23Feb2019)

# Dose Analysis – RCT

| 5-year Mortality Rate |              |              |              |              |
|-----------------------|--------------|--------------|--------------|--------------|
| Dose Group 1          | Dose Group 2 | Dose Group 3 | Dose Group 4 | Dose Group 5 |
| 11.5%                 | 13.6%        | 13.4%        | 20.0%        | 13.2%        |
| p=0.72                |              |              |              |              |



**No impact of Zilver PTX paclitaxel dose on mortality rate**


# Causes of Death Through 5 Years – RCT and BMS

| Cause            | RCT – PTX<br>(n=336) | RCT – PTA / BMS<br>(n=143) | p-value | Zilver BMS Study*<br>(n=110) |
|------------------|----------------------|----------------------------|---------|------------------------------|
| Cardiovascular   | 4.8%                 | 5.6%                       | 0.66    | 4.5%                         |
| Cancer           | 4.8%                 | 1.4%                       | 0.11    | 6.4%                         |
| Pulmonary        | 1.8%                 | 1.4%                       | > 0.99  | 1.8%                         |
| Stroke           | 0.6%                 | 0.7%                       | > 0.99  | 0.0%                         |
| Trauma           | 0.0%                 | 1.4%                       | 0.09    | 0.0%                         |
| GI               | 0.3%                 | 0.0%                       | > 0.99  | 0.9%                         |
| Multiple/Unknown | 2.1%                 | 1.4%                       | > 0.99  | 0.9%                         |

**No increased rate of cardiovascular, cancer, or other cause  
of death for Zilver PTX compared to PTA or BMS**




# Japan Post-Market Studies – Zilver PTX and BMS



- No exclusion criteria
  - Challenging patient population, including CLI patients
- 904 Zilver PTX patients
  - 5-year follow-up
- 190 BMS patients
  - 3-year follow-up
  - Separate study, not randomized
- No significant difference in mortality ( $p=0.92$ )
- Same mortality rate of 5.1% per year for PTX & BMS
  - Linear from 0-3 and 3-5 years


# Japan Post-Market Studies – Zilver PTX and BMS Claudicants



| Mortality |       |      |
|-----------|-------|------|
|           | PTX   | BMS  |
| 3 years   | 10.0% | 8.8% |
| 5 years   | 19.7% | -    |

No significant difference  
in mortality ( $p=0.57$ )

# Japan Post-Market Studies – Zilver PTX and BMS CLI Patients



| Mortality |       |       |
|-----------|-------|-------|
|           | PTX   | BMS   |
| 3 years   | 37.7% | 32.2% |
| 5 years   | 51.3% | -     |

No significant difference  
in mortality ( $p=0.72$ )

# Covariate Analysis – Japan

- Cox proportional hazards model to evaluate covariates
  - No significant difference between Zilver PTX and BMS ( $p=0.39$ )

| Cox Model Results    |                      |
|----------------------|----------------------|
| Covariate            | Multivariate p-value |
| Age                  | <0.0001              |
| Claudication/CLI     | <0.0001              |
| Hypercholesterolemia | 0.0005               |
| Gender               | 0.003                |
| Diabetes             | 0.04                 |
| Carotid disease      | 0.06                 |
| PTX vs. BMS          | 0.39                 |
| Smoking              | 0.45                 |
| Hypertension         | 0.46                 |
| Lesion length        | 0.80                 |
| Pulmonary disease    | 0.90                 |

# Dose Analysis – Japan

| 5-year Mortality Rate |              |              |              |              |
|-----------------------|--------------|--------------|--------------|--------------|
| Dose Group 1          | Dose Group 2 | Dose Group 3 | Dose Group 4 | Dose Group 5 |
| 17.4%                 | 23.9%        | 16.1%        | 21.3%        | 21.5%        |
| p=0.41                |              |              |              |              |



**No impact of Zilver PTX paclitaxel dose on mortality rate**

# Causes of Death Through 5 Years – RCT & Japan

| Cause            | RCT – PTX<br>(n=336) | RCT – PTA / BMS<br>(n=143) | Japan – PTX<br>(n=904)* |
|------------------|----------------------|----------------------------|-------------------------|
| Cardiovascular   | 4.8%                 | 5.6%                       | 6.1%                    |
| Cancer           | 4.8%                 | 1.4%                       | 2.9%                    |
| Pulmonary        | 1.8%                 | 1.4%                       | 2.7%                    |
| Stroke           | 0.6%                 | 0.7%                       | 1.5%                    |
| Trauma/Accident  | 0.0%                 | 1.4%                       | 0.2%                    |
| GI               | 0.3%                 | 0%                         | 0.2%                    |
| Infection        | 0%                   | 0%                         | 0.2%                    |
| Renal            | 0%                   | 0%                         | 0.8%                    |
| Multiple/Unknown | 2.1%                 | 1.4%                       | 5.9%                    |

# Mortality Rates from Literature

## 3-YEAR MORTALITY

Cook  
PTX

Cook  
BMS


PAD Mortality  
Literature Review

## 5-YEAR MORTALITY

Cook  
PTX


Cook  
BMS

PAD Mortality  
Literature Review



# No Increased Long-term Mortality with DES

- 51,456 patients
  - 47,351 BMS
  - 4,105 DES (Zilver PTX)
- Similar mortality for BMS and DES through 4.1 years
  - Overall adjusted  $p=0.53$
  - Without CLI adjusted  $p=0.95$
  - With CLI adjusted  $p=0.32$



Secemsky E, et al. J Am Coll Cardiol. E-pub ahead of print 01March2019. doi <https://doi.org/10.1016/j.jacc.2019.02.020>



# Conclusions

- Conclusion of Katsanos K, et al. was not based on patient-level data
- Patient-level analysis of RCT data shows no increased long-term mortality risk with Zilver PTX compared to PTA and BMS
  - Covariate analysis supports no significant difference
  - No impact of Zilver PTX paclitaxel dose on mortality rate
  - No significant differences in causes of death
- Mortality rates for the Zilver PTX stent are consistent with rates reported in literature for PAD patients
- Japan data confirm RCT findings showing no increased long-term mortality risk with Zilver PTX compared to BMS
- Cook will continue to work with global regulatory authorities and independent physician led groups to evaluate safety using patient-level data